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
Fig.2. Cytogram of Total SWS vs total Yellow Green fluorescence (gated cluster indicate heterotrophic nanoflagellates)

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1. Dubelaar, G. B. J., A. Groenewegen, et al. (1989). "Optical plankton analyser: a flow cytometer for plankton analysis, II: Specifications." *Cytometry* 10(5): 529-539.
2. Peeters, J. C. H., G. B. J. Dubelaar, et al. (1989). "Optical plankton analyser: A flow cytometer for plankton analysis, I: Design considerations." *Cytometry* 10(5): 522-528.
3. Jonker, R. R., G. B. J. Dubelaar, et al. (1994). "The European Optical Plankton Analyser: A high dynamic range flow cytometer." *Scientia Marina*.
4. Jonker, R. R., J. T. Meulemans, et al. (1995). "Flow cytometry: A powerful tool in analysis of biomass distributions in phytoplankton " *Water Sci. Technol.* 32(4): 177-182.
5. Dubelaar, G. B. J., P. L. Gerritzen, et al. (1999). "Design and first results of CytoBuoy: A wireless flow cytometer for in situ analysis of marine and fresh waters." *Cytometry* 37(4): 247-254.
6. Wilkins, M. F., L. Boddy, et al. (1999). "Identification of Phytoplankton from Flow Cytometry Data by Using Radial Basis Function Neural Networks." *Appl. Environ. Microbiol.* 65(10): 4404-4410.
7. Dubelaar, G. B. J. and P. L. Gerritzen (2000). "CytoBuoy: a step forward towards using flow cytometry in operational oceanography." *Sci. Mar. (Barc.)* 64(2): 255-265.
8. Dubelaar, G. B. J. and R. R. Jonker (2000). "Flow cytometry as a tool for the study of phytoplankton." *Scientia Marina* 64(2): 135-156
9. Woodd-Walker, R. S., C. P. Gallienne, et al. (2000). "A test model for optical plankton counter (OPC) coincidence and a comparison of OPC-derived and conventional measures of plankton abundance." *J. Plankton Res.* 22(3): 473-483.
10. Reynolds, C. S., V. Huszar, et al. (2002). "Towards a functional classification of the freshwater phytoplankton." *J. Plankton Res.* 24(5): 417-428.
11. Cunninghama, A., D. McKee, et al. (2003). "Fine-scale variability in phytoplankton community structure and inherent optical properties measured from an autonomous underwater vehicle " *J. Mar. Syst.* 43(1-2): 51-59.
12. Dubelaar, G. B. J. and P. J. F. Geerders (2004). "Innovative technologies to monitor plankton dynamics " *Sea Technol.* 45(8): 15-21.
13. Dubelaar, G. B. J., P. J. F. Geerders, et al. (2004). "High frequency monitoring reveals phytoplankton dynamics." *J. Environ. Monit.* 6: 946-952.
14. Takabayashi, M., F. P. Wilkerson, et al. (2005). "Response Of Glutamine Synthetase Gene Transcription And Enzyme Activity To External Nitrogen Sources In The Diatom *Skeletonema Costatum* (Bacillariophyceae)." *J. Phycol.* 41(1): 84-94.
15. Takabayashi, M., K. Lew, et al. (2006). "The effect of nutrient availability and temperature on chain length of the diatom, *Skeletonema costatum*." *J. Plankton Res.* 28(9): 831-840.
16. Dubelaar, G. B. J., R. Casotti, et al. (2007). "Phytoplankton and their analysis by flow cytometry." *Flow Cytometry with Plant Cells*: 287-322.

17. 17. Thyssen, M., D. Mathieu, et al. (2008). "Short-term variation of phytoplankton assemblages in Mediterranean coastal waters recorded with an automated submerged flow cytometer." *J. Plankton Res.* 30(9): 1027-1040.
18. 18. Thyssen, M. T., N. Garcia, et al. (2008). "Sub meso scale phytoplankton distribution in the north east Atlantic surface waters determined with an automated flow cytometer." *Biogeosciences Discuss* 5: 2471-2503.
19. 19. Thyssen, M. T., G. A. Zubkov, et al. (2008). "The emergence of automated high-frequency flow cytometry: revealing temporal and spatial phytoplankton variability." *J. Plankton Res.* 30(3): 333-343.
20. 20. Barofsky, A., P. Simonelli, et al. (2009). "Growth phase of the diatom *Skeletonema marinoi* influences the metabolic profile of the cells and the selective feeding of the copepod *Calanus* spp." *J. Plankton Res.* 32(3): 263-272.
21. 21. Pereira, G. C., A. Granato, et al. (2009). "Viriplankton Abundance in Trophic Gradients of an Upwelling Field." *Brazilian Journal of Microbiology* 40: 857-865.
22. 22. Van Donk, E., S. Cerbin, et al. (2009). "The effect of a mixotrophic chrysophyte on toxic and colony-forming cyanobacteria." *Freshwater Biology* 54(9): 1843-1855.
23. 23. Barofsky, A., P. Simonelli, et al. (2010). "Growth phase of the diatom *Skeletonema marinoi* influences the metabolic profile of the cells and the selective feeding of the copepod *Calanus* spp." *J. Plankton Res.* 32(3): 263-272.
24. 24. Hansen, B. W., H. H. Jakobsen, et al. (2010). "Swimming behavior and prey retention of the polychaete larvae *Polydora ciliata*." *Journal of Experimental Biology*(213): 3237-3246.
25. 25. Pereira, G., A. de Figueiredo, et al. (2010). "Assessing the ecological status of plankton in Anjos Bay: a flow cytometry approach." *Biogeosciences Discussions* 7: 6243-6264.
26. 26. Thyssen, M., B. Beker, et al. (2011). "Phytoplankton distribution during two contrasted summers in a Mediterranean harbour: combining automated submersible flow cytometry with conventional techniques." *Environmental Monitoring and Assessment* 173(1): 1-16.
27. 27. Thyssen, M. and M. Denis (2011). *Temporal and Spatial High-Frequency Monitoring of Phytoplankton by Automated Flow Cytometry and Pulse-Shape Analysis*. *Global Change: Mankind-Marine Environment Interactions*. H.-J. Ceccaldi, I. Dekeyser, M. Girault and G. Stora, Springer Netherlands: 293-298.
28. 28. Vidoudez, C., J. C. Nejstgaard, et al. (2011). "Dynamics of Dissolved and Particulate Polyunsaturated Aldehydes in Mesocosms Inoculated with Different Densities of the Diatom *Skeletonema marinoi*." *Marine Drugs* 9(3): 345-358.
29. 29. Malkassian, A., D. Nerini, et al. (2011). "Functional analysis and classification of phytoplankton based on data from an automated flow cytometer." *Cytometry Part A* 79A(4): 263-275.
30. 30. G.C. Pereira, N.F.F. Ebecken *Combining In Situ Flow Cytometry And Artificial Neural Networks For Aquatic Systems Monitoring Expert Systems with Applications*
31. 31. Juliano C. Cury, Fabio V. Araujo, Sergio A. Coelho-Souza, Raquel S. Peixoto, Joana A. L. Oliveira, Henrique F. Santos, Alberto M. R. Dávila, and Alexandre S. Rosado *Microbial Diversity Of A Brazilian Coastal Region Influenced By An Upwelling System And Anthropogenic Activity* *LoS One*. 2011; 6(1): e16553. Published online 2011 January 27. doi: 10.1371/journal.pone.0016553.

32. 32. Thyssen, M., Ferreyra, G., Moreau, S., Schloss S., Denis, M., Demers, S. Combined Effects Of Ultraviolet Radiations B And Temperature Increase On Phytoplankton Dynamics And Cell Cycle Using Pulse Shape Recording Flow Cytometry. *Journal of Experimental Marine Biology and Ecology*, under acceptance.
33. 33. Francesco Pomati, Jukka Jokela, Marco Simona, Mauro Veronesi, and B. W. Ibelings An Automated Platform for Phytoplankton Ecology and Aquatic Ecosystem Monitoring *Environ. Sci. Technol.*, 2011, 45 (22), pp 9658–9665
34. 34. Q. Zhou, W. Chen, H. Zhang et al. A Flow Cytometer Based Protocol for Quantitative Analysis of Bloom-forming Cyanobacteria (*Microcystis*) in Lake Sediments *Journal of Environmental Sciences* 2012, 24(9) 1709–1716
35. 35. Johanna Bergkvist, Peter Thor, Hans Henrik Jakobsen, Sten-Åke Wängberg and Erik Selander Grazer-induced chain length plasticity reduces grazing risk in a marine diatom *Limnol. Oceanogr.*, 57(1), 2012, 318–324
36. 36. Andrea Gerech, Ylenia Carotenuto, Adrianna Ianora, Giovanna Romano, Angelo Fontana, Giuliana d'Ippolito, Hans H. Jakobsen, Jens C. Nejstgaard Oxylin production during a mesocosm bloom of *Skeletonema marinoi* *Journal of Experimental Marine Biology and Ecology*, August 2013, Volume 446, 159-165
37. 37. Francesco Pomati, Nathan J. B. Kraft, Thomas Posch, Bettina Eugster, Jukka Jokela, Bas W. Ibelings Individual Cell Based Traits Obtained by Scanning Flow-Cytometry Show Selection by Biotic and Abiotic Environmental Factors during a Phytoplankton Spring Bloom *PLOS ONE*, August 2013, Volume 8, Issue 8, e71677
38. M. Thyssen, S. Alvain, a. Lefèbvre, D. Dessailly, M. Rijkeboer, N. Guiselin, V. Creach, L.-F. Artigas, Phytoplankton community structure in the North Sea: coupling between remote sensing and automated in situ analysis at the single cell level, *Biogeosciences Discuss.* 11 (2014) 15621–15662.
39. M. Latasa, J. Henjes, R. Scharek, P. Assmy, R. Röttgers, V. Smetacek, Progressive decoupling between phytoplankton growth and microzooplankton grazing during an iron-induced phytoplankton bloom in the Southern Ocean (EIFEX), *Mar. Ecol. Prog. Ser.* 513 (2014) 39–50.
40. M. Dugenne, M. Thyssen, D. Nerini, C. Mante, J.C. Poggiale, N. Garcia, F. Garcia, G.J. Grégori, Consequence of a sudden wind event on the dynamics of a coastal phytoplankton community: An insight into specific population growth rates using a single cell high frequency approach, *Front. Microbiol.* 5 (2014) 1–14.
41. M. Thyssen, G.J. Grégori, J.-M. Grisoni, M.L. Pedrotti, L. Mousseau, L.F. Artigas, S. Marro, N. Garcia, O. Passafiume, M.J. Denis, Onset of the spring bloom in the northwestern Mediterranean Sea: influence of environmental pulse events on the in situ hourly-scale dynamics of the phytoplankton community structure, *Front. Microbiol.* 5 (2014) 387.
42. S. Fontana, J. Jokela, F. Pomati, Opportunities and challenges in deriving phytoplankton diversity measures from individual trait-based data obtained by scanning flow-cytometry, *Front. Microbiol.* 5 (2014) 1–12.
43. E.D. Cooper, B. Bentlage, T.R. Gibbons, T.R. Bachvaroff, C.F. Delwiche, Metatranscriptome profiling of a harmful algal bloom, *Harmful Algae.* 37 (2014) 75–83.

44. G.C. Pereira, L.P. Andrade, R.P. Espíndola, N.F.F. Ebecken, Structural Analysis and Static Simulation of Coastal Planktonic Networks, *JILSA* (2014) 113–124.
45. M. Denis, Recent advances in assessing the Dynamics of phytoplankton assemblages by high frequency analysis at the single cell level, *CIBS Meeting Proceedings*, (2014) 77–85.
46. S. Glüge, F. Pomati, C. Albert, P. Kauf, T. Ott, The challenge of clustering flow cytometry data from phytoplankton in lakes, In *Nonlinear Dynamics of Electronic Systems*, Springer International Publishing, (2014) 379-389.
47. M.N. McFarland, Impact of Phytoplankton Community Structure and Function on Marine Particulate Optical Properties, *Open Access Dissertation* (2014).
48. K.R. Owen, Flow Cytometric Investigation of the Size Spectrum of North Sea Phytoplankton Communities, *Open Access Dissertation* (2014).
49. L. Duforêt-Gaurier, W. Moutier, N. Guiselin, M. Thyssen, G. Dubelaar, X. Mériaux, L. Courcot, D. Dessailly, H. Loisel, Determination of backscattering cross section of individual particles from cytometric measurements: a new methodology, *Opt. Express*. 23 (2015) 31510-31533.
50. N.F.F. Ebecken, G.C. Pereira, Ecological Networks : From Structural Analysis to Synchronization, 9 (2015) 938–942.
51. S. Fontana, O.L. Petchey, F. Pomati, Individual-level trait diversity concepts and indices to comprehensively describe community change in multidimensional trait space, *Funct. Ecol.* 30 (2015) 808–818.
52. M. Thyssen, S. Alvain, A. Lefèbvre, D. Dessailly, M. Rijkeboer, N. Guiselin, V. Creach, L.-F. Artigas, High-resolution analysis of a North Sea phytoplankton community structure based on in situ flow cytometry observations and potential implication for remote sensing, *Biogeosciences*. 12 (2015) 4051–4066.
53. M.N. McFarland, J. Rines, J. Sullivan, P. Donaghay, Impact of phytoplankton size and physiology on particulate optical properties determined with scanning flow cytometry, *Mar. Ecol. Prog. Ser.* 531 (2015) 43–61.
54. S.A. Coelho-Souza, F. V. Araujo, J.C. Cury, H.E. Jesus, G.C. Pereira, J.R.D. Guimaraes, R.S. Peixoto, A.M.R. Davila, A.S. Rosado, Bacterial and archaeal communities variability associated with upwelling and anthropogenic pressures in the protection area of Arraial do Cabo (Cabo Frio region - RJ), *An. Acad. Bras. Cienc.* 87 (2015) 1737–1750.
55. S. Bonato, U. Christaki, A. Lefebvre, F. Lizon, M. Thyssen, L.F. Artigas, High spatial variability of phytoplankton assessed by flow cytometry, in a dynamic productive coastal area, in spring: The eastern English Channel, *Estuar. Coast. Shelf Sci.* 154 (2015) 214–223.
56. M. Dugenne, M. Thyssen, N. Garcia, N. Mayot, G. Bernard, Monitoring of a Potential Harmful Algal Species in the Berre Lagoon by Automated In Situ Flow Cytometry, In *Marine Productivity: Perturbations and Resilience of Socio-ecosystems*, Springer International Publishing, (2015) 117–127.
57. A. Golabgir, D. Ehgartner, L. Neutsch, A.E. Posch, P. Sagmeister, C. Herwig, Imaging Flow Cytometry and High-Throughput Microscopy for Automated Macroscopic Morphological Analysis of Filamentous Fungi, *InGenetic Transformations Systems in Fumgi*, Volume 2, Springer International Publishing, (2015) 201–210.
58. D. Ehgartner, J. Fricke, A. Schröder, C. Herwig, At-line determining spore germination of *Penicillium chrysogenum* bioprocesses in complex media, *Appl. Microbiol. Biotechnol.* 100 (2016) 8923–8930.

59. D. Ehgartner, C. Herwig, L. Neutsch, At-line determination of spore inoculum quality in *Penicillium chrysogenum* bioprocesses, *Appl. Microbiol. Biotechnol.* 100 (2016) 5363–5373.
60. S. Bonato, E. Breton, M. Didry, F. Lizon, V. Cornille, E. Lécuyer, U. Christaki, L.F. Artigas, Spatio-temporal patterns in phytoplankton assemblages in inshore-offshore gradients using flow cytometry: A case study in the eastern English Channel, *J. Mar. Syst.* 156 (2016) 76–85.
61. S. Fontana, M. K. Thomas, M. Moldoveanu, P. Spaak, F. Pomati, Individual-level trait diversity predicts phytoplankton community properties better than species richness or evenness, *The ISME Journal* (2017) 1-11
62. T. Silovic, G. Grégori, M. Dugenne, M. Thyssen, F. Calendreau, T. Cossart, H. Kools, G. Dubelaar, M. Denis, A new automated flow cytometer for high frequency in situ characterisation of heterotrophic microorganisms and their dynamics in aquatic ecosystems, IMEKO International Conference on Metrology for The Sea, Naples, Italy, October 11-13, 2017
63. D. Ehgartner, C. Herwig, J. Fricke, Morphological analysis of the filamentous fungus *Penicillium chrysogenum* using flow cytometry—the fast alternative to microscopic image analysis, *Appl. Microbiol. Biotechnol.* 101 (2017) 7675-7688
64. G.C. Pereira, A.R. Figueiredo, N.F.F. Ebecken, Using in situ flow cytometry images of ciliates and dinoflagellates for aquatic system monitoring, *Braz. J. Biol.* (2017)
65. L. Haraguchi, H. H. Jakobsen, N. Lundholm, J. Carstensen, Monitoring natural phytoplankton communities: a comparison between traditional methods and pulse-shape recording flow cytometry, *Aquat. Microb. Ecol.* 80 (2017) 77-92
66. W. Moutier, L. Duforêt-Gaurier, M. Thyssen, H. Loisel, X. Mériaux, L. Courcot, D. Dessailly, A.-H. Rêve, G. Grégori, S. Alvain, A. Barani, L. Brutier, M. Dugenne, Evolution of the scattering properties of phytoplankton cells from flow cytometry measurements, *PLoS ONE* (2017)
67. E. Breton, U. Christaki, S. Bonato, M. Didry, L.F. Artigas, Functional trait variation and nitrogen use efficiency in temperate coastal phytoplankton, *Mar. Ecol. Prog. Ser.* 563 (2017) 35-49
68. H. Tan, T. Oishi, A. Tanaka, R. Doerffer, Y. Tan, Chlorophyll-a specific volume scattering function of phytoplankton, *Optics EXPRESS* 25 (2017) A564-A573
69. F. Pomati, J. Jokela, S. Castiglioni, M. K. Thomas, L. Nizzetto, Water-borne pharmaceuticals reduce phenotypic diversity and response capacity of natural phytoplankton communities, *PLoS ONE* (2017)
70. J. Taucher, M. Haunost, T. Boxhammer, L. T. Bach, M. Alguero-Muñiz, U. Riebesell, Influence of ocean acidification on plankton community structure during a winter-to-summer succession: An imaging approach indicates that copepods can benefit from elevated CO₂ via indirect food web effects, *PLoS ONE* (2017)
71. J. Taucher, J. Arístegui, L.T. Bach, W. Guan, M.F. Montero, A. Nauendor, E.P. Achterberg, U. Riebesell, Response of Subtropical Phytoplankton Communities to Ocean Acidification Under Oligotrophic Conditions and During Nutrient Fertilization, *Front. Mar. Sci.* (doi: 10.3389/fmars.2018.00330)
72. M.K. Thomas, S. Fontana, M. Reyes, F. Pomati, Quantifying cell densities and biovolumes of phytoplankton communities and functional groups using scanning flow cytometry, machine learning and unsupervised clustering, *PLoS ONE* (2018)
73. M.K. Thomas, S. Fontana, M. Reyes, M. Kehoe, F. Pomati, The predictability of a lake phytoplankton community, over time-scales of hours to years, *Ecol. Lett.* (2018)

74. K. Leblanc, B. Quéguiner, F. Diaz, V. Cornet, M. Michel-Rodriguez, X. Durrieu de Madron, C. Bowler, S. Malviya, M. Thyssen, G. Grégori, M. Rembauville, O. Grosso, J. Poulain, C. de Vargas, M. Pujo-Pay, P. Conan, Nanoplanktonic diatoms are globally overlooked but play a role in spring blooms and carbon export, *Nature Comm.* 9 (2018)
75. J. Agagliate, R. Röttgers, M. S. Twardowski, D. McKee, Evaluation of a flow cytometry method to determine size and real refractive index distributions in natural marine particle populations, *Appl. Optics* 57 (2018) 1705-1716
76. J. Agagliate, I. Lefering, D. McKee, Forward modeling of inherent optical properties from flow cytometry estimates of particle size and refractive index, *Appl. Optics* 57 (2018) 1777-1788
77. L. T. Bach, K. T. Lohbeck, T. B. H. Reusch, U. Riebesell, Rapid evolution of highly variable competitive abilities in a key phytoplankton species, *Nature Ecol. Evol.* (2018)

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